


| | | Mini | | | | Regular | | | | | | | | | | | | | | | | | |
|--------------|--------|---------------------|---------------|---------------------|--------------|-------------------------|--------------|---------------------|--------------|-------------------------|--------------|-------------------------|--------------|-------------------------|--------------|-------------------------|--------------|-------------------------|--------------|---------------------|--------------|-------------------------|--------|
| | | Φ3.5 | | Φ4.0 | | Φ4.5 | | Φ5.0 | | Φ6.0 | | Φ7.0 | | | | | | | | | | | |
| 对应编码 | | TS3M35 | | | | TS3S40 | | | | TS3S45 | | | | TS3S50 | | TS3S60 | | | | TS3S70 | | | |
| S-L-A 种植体 | L | 7mm | 311010 | | | | 312010 | | | | 313010 | | | | 314020 | | 315020 | | | | 316020 | | 7mm |
| | | 8.5mm | 311020 | | | | 312020 | | | | 313020 | | | | 314030 | | 315030 | | | | 316030 | | 8.5mm |
| | | 10mm | 311030 | | | | 312030 | | | | 313030 | | | | 314040 | | 315040 | | | | 316040 | | 10mm |
| | | 11.5mm | 311040 | | | | 312040 | | | | 313040 | | | | 314050 | | 315050 | | | | 316050 | | 11.5mm |
| | | 13mm | 311050 | | | | 312050 | | | | 313050 | | | | 314060 | | 315060 | | | | 316060 | | 13mm |
| 15mm | | 311050 | | | | 312060 | | | | 313060 | | | | 314070 | | 315060 | | | | 316060 | | 15mm | |
| 三棱钻 | | | | | | | | | | | | | | | | | | | | | | | |
| | | Φ2.0 351030 | | Φ2.0 351030 | | Φ2.0 351030 | | Φ2.0 351030 | | Φ2.0 351030 | | | | | | | | | | | | | |
| 先锋钻 | 种植体 L | | | | | | | | | | | | | | | | | | | | | | |
| | | Φ2.0 | | Φ2.0 | | Φ2.0 | | Φ2.0 | | Φ2.0 | | Φ2.2 352500 | | Φ2.2 352500 | | Φ2.2 352500 | | Φ2.2 352500 | | | | | |
| | | 7mm | 351080 | | | | 351080 | | | | 351080 | | | | 7mm | | 351080 | | | | 7mm | | |
| | | 8.5mm | 351110 | | | | 351110 | | | | 351110 | | | | 8.5mm | | 351110 | | | | 8.5mm | | |
| | | 10mm | 351140 | | | | 351140 | | | | 351140 | | | | 10mm | | 351140 | | | | 10mm | | |
| | | 11.5mm | 351170 | | | | 351170 | | | | 351170 | | | | 11.5mm | | 351170 | | | | 11.5mm | | |
| 13mm | | 351250 | | | | 351250 | | | | 351250 | | | | 13mm | | 351250 | | | | 13mm | | | |
| 15mm | | 351380 | | | | 351380 | | | | 351380 | | | | 15mm | | 351380 | | | | 15mm | | | |
| 扩孔钻 | 种植体 L | | | | | | | | | | | | | | | | | | | | | | |
| | | Φ3.0(松质骨) | | F3.5(标准骨) | | F3.5(松质骨) | | F4.0(标准骨) | | F3.5 | | F4.5(松质骨) | | F6.0(松质骨) (备注：扩孔至刻度线下缘) | | F6.0 (备注：扩孔至刻度线上缘) | | F7.0(松质骨) (备注：扩孔至刻度线下缘) | | | | | |
| | | 7mm | 351180/351300 | | | | 351700 | | 351760 | | 351700 | | 351820 | | 7mm | | 352530 | | | | 7mm | | |
| | | 8.5mm | 351180/351300 | | | | 351710 | | 351760 | | 351710 | | 351830 | | 8.5mm | | 352540 | | | | 8.5mm | | |
| | | 10mm | 351180/351300 | | | | 351720 | | 351780 | | 351720 | | 351840 | | 10mm | | 352550 | | | | 10mm | | |
| | 11.5mm | 351180/351300 | | | | 351730 | | 351790 | | 351730 | | 351850 | | 11.5mm | | 352560 | | | | 11.5mm | | | |
| | 13mm | 351180/351300 | | | | 351740 | | 351800 | | 351740 | | 351860 | | 13mm | | 352570 | | | | 13mm | | | |
| | 15mm | 351180/351300 | | | | 351750 | | 351810 | | 351750 | | 351870 | | 15mm | | 352570 | | | | 15mm | | | |
| | 种植体 L | | | | | | | | | | | | | | | | | | | | | | |
| | | F4.5(标准骨) | | F5.0(标准骨) | | F6.0(标准骨) (备注：扩孔至刻度线上缘) | | F6.0 (备注：扩孔至刻度线上缘) | | F7.0(标准骨) (备注：扩孔至刻度线上缘) | | F6.0(标准骨) (备注：扩孔至刻度线上缘) | | F6.0 (备注：扩孔至刻度线上缘) | | F7.0(标准骨) (备注：扩孔至刻度线上缘) | | F6.0(标准骨) (备注：扩孔至刻度线上缘) | | F6.0 (备注：扩孔至刻度线上缘) | | F7.0(标准骨) (备注：扩孔至刻度线上缘) | |
| 7mm | | | | | | 351820 | | 351880 | | 351820 | | 7mm | | 352530 | | | | 352530 | | 352590 | | 7mm | |
| 8.5mm | | | | | | 351830 | | 351890 | | 351830 | | 8.5mm | | 352540 | | | | 352540 | | 352600 | | 8.5mm | |
| 10mm | | | | | | 351840 | | 351900 | | 351840 | | 10mm | | 352550 | | | | 352550 | | 352610 | | 10mm | |
| 11.5mm | | | | | 351850 | | 351910 | | 351850 | | 11.5mm | | 352560 | | | | 352560 | | 352620 | | 11.5mm | | |
| 13mm | | | | | 351860 | | 351920 | | 351860 | | 13mm | | 352570 | | | | 352570 | | 352630 | | 13mm | | |
| 15mm | | | | | 351870 | | 351930 | | 351870 | | 15mm | | 352570 | | | | 352570 | | 352630 | | 15mm | | |
| 皮质骨钻 | | | | | | | | | | | | | | | | | | | | | | | |
| | | F3.5(硬质骨) 351610 | | F4.0(硬质骨) 351620 | | F4.5(硬质骨) 351630 | | F5.0(硬质骨) 351640 | | F6.0(硬质骨) 352650 | | F7.0(硬质骨) 352660 | | F6.0(硬质骨) 352650 | | F7.0(硬质骨) 352660 | | F6.0(硬质骨) 352650 | | F7.0(硬质骨) 352660 | | | |
| 封闭帽 愈合基台 | 封闭螺丝 | Φ4.0 | | Φ4.5 | | Φ4.0 | | Φ5.0 | | Φ5.5 | | Φ6.0 | | Φ7.0 | | Φ6.0 | | Φ7.0 | | Φ7.0 | | | |
| | | 321010(H0.4) | 323010(H3.0) | 323110(H3.0) | 323310(H3.0) | 322010(H0.4) | 324010(H3.0) | 324110(H3.0) | 324510(H3.0) | 324210(H3.0) | 324310(H3.0) | 321020(H2.0) | 323020(H4.0) | 323120(H4.0) | 323340(H4.0) | 322020(H2.0) | 324020(H4.0) | 324120(H4.0) | 324520(H4.0) | 324220(H4.0) | 324320(H4.0) | | |
| | | 323030(H5.0) | 323130(H5.0) | 323320(H5.0) | 323330(H5.0) | 324030(H5.0) | 324130(H5.0) | 324530(H5.0) | 324230(H5.0) | 324330(H5.0) | 324340(H5.0) | 323040(H7.0) | 323140(H7.0) | 323330(H7.0) | 323330(H7.0) | 324040(H7.0) | 324140(H7.0) | 324540(H7.0) | 324240(H7.0) | 324340(H7.0) | 324340(H7.0) | | |
| | | 323040(H7.0) | 323140(H7.0) | 323330(H7.0) | 323330(H7.0) | 324040(H7.0) | 324140(H7.0) | 324540(H7.0) | 324240(H7.0) | 324340(H7.0) | 324340(H7.0) | 323040(H7.0) | 323140(H7.0) | 323330(H7.0) | 323330(H7.0) | 324040(H7.0) | 324140(H7.0) | 324540(H7.0) | 324240(H7.0) | 324340(H7.0) | 324340(H7.0) | | |
| | | 323040(H7.0) | 323140(H7.0) | 323330(H7.0) | 323330(H7.0) | 324040(H7.0) | 324140(H7.0) | 324540(H7.0) | 324240(H7.0) | 324340(H7.0) | 324340(H7.0) | 323040(H7.0) | 323140(H7.0) | 323330(H7.0) | 323330(H7.0) | 324040(H7.0) | 324140(H7.0) | 324540(H7.0) | 324240(H7.0) | 324340(H7.0) | 324340(H7.0) | | |

Mini

Regular

临时修复体/保护帽



临时基台


| H | GH | Type | D | Φ4.0 |
|----|-----|------|---|--------|
| 10 | 1.0 | 六角 | | 371010 |
| | | 非六角 | | 371020 |
| | 3.0 | 六角 | | 371030 |
| | | 非六角 | | 371040 |

* 建议螺栓拧紧扭矩: 20-30Ncm

多功能基台保护帽



377040



临时基台

| H | GH | Type | D | Φ4.5 |
|----|-----|------|---|--------|
| 10 | 1.0 | 六角 | | 371050 |
| | | 非六角 | | 371060 |
| | 3.0 | 六角 | | 371070 |
| | | 非六角 | | 371080 |

* 建议螺栓拧紧扭矩: 20-35Ncm

多功能基台保护帽







377040

愈合帽



325010 325020 325030

印模转移



368010

印模帽

代体



361010



362140



362010



362140



368020

修复方式

粘固固位
螺栓粘固混合固位

粘固固位
螺栓粘固混合固位

螺栓固位

粘固固位
螺栓粘固混合固位

粘固固位
螺栓粘固混合固位

螺栓固位


球帽固位

基台



| H | GH | D | Φ4.5 |
|-----|-----|----|--------|
| 5.5 | 1.0 | 六角 | 331010 |
| | | | 331020 |
| | | | 331030 |
| | | | 331040 |
| | | | 331050 |
| 7.0 | 1.0 | | 331510 |
| | | | 331520 |
| | | | 331530 |
| | | | 331540 |
| | | | 331550 |

* 建议螺栓拧紧扭矩: 30Ncm




| H | GH | Type | D | Φ4.5 |
|-----|-----|------|---|--------|
| 8.0 | 2.0 | A | | 333010 |
| | | B | | 333020 |
| | 4.0 | A | | 333040 |
| | | B | | 333050 |

* 建议螺栓拧紧扭矩: 30Ncm




| GH | D | Φ4.8 |
|-----|---|--------|
| 1.0 | | 337010 |
| 2.0 | | 337020 |
| 3.0 | | 337030 |
| 4.0 | | 337040 |
| 5.0 | | 337050 |

* 建议螺栓拧紧扭矩: 30Ncm



| A | GH | D | Φ4.8 |
|-----|-----|---|--------|
| 17° | 2.5 | | 337130 |
| | 3.0 | | 337140 |
| | 4.0 | | 337150 |
| 30° | 3.5 | | 337160 |
| | 4.0 | | 337170 |
| | 5.0 | | 337180 |

* 建议螺栓拧紧扭矩: 25Ncm



| H | GH | D | Φ4.5 | Φ5.0 | Φ6.0 | Φ7.0 |
|-----|-----|----|--------|--------|--------|--------|
| 5.5 | 1.0 | 六角 | 331060 | 331160 | 331310 | 331410 |
| | | | 331070 | 331170 | 331320 | 331420 |
| | | | 331080 | 331180 | 331330 | 331430 |
| | | | 331090 | 331190 | 331340 | 331440 |
| | | | 331100 | 331200 | 331350 | 331450 |
| 7.0 | 1.0 | | 331110 | 331210 | 331360 | 331460 |
| | | | 331120 | 331220 | 331370 | 331470 |
| | | | 331130 | 331230 | 331380 | 331480 |
| | | | 331140 | 331240 | 331390 | 331490 |
| | | | 331150 | 331250 | 331400 | 331500 |

* 建议螺栓拧紧扭矩: 35Ncm




| H | GH | Type | Φ5.0 | Φ6.0 |
|-----|-----|------|--------|--------|
| 8.0 | 2.0 | A | 333070 | 333130 |
| | | B | 333080 | 333140 |
| | 4.0 | A | 333100 | 333160 |
| | | B | 333110 | 333170 |

* 建议螺栓拧紧扭矩: 35Ncm



| GH | D | Φ4.8 |
|-----|---|--------|
| 1.0 | | 337070 |
| 2.0 | | 337080 |
| 3.0 | | 337090 |
| 4.0 | | 337100 |
| 5.0 | | 337110 |

* 建议螺栓拧紧扭矩: 30Ncm



| A | GH | D | Φ4.8 |
|-----|-----|---|--------|
| 17° | 2.5 | | 337210 |
| | 3.0 | | 337200 |
| | 4.0 | | 337210 |
| | 3.5 | | 337220 |
| 30° | 4.0 | | 337230 |
| | 5.0 | | 337240 |

* 建议螺栓拧紧扭矩: 30Ncm



| 穿眼高度 | 迷你 | | 标准 |
|------|--------|--------|--------|
| | 1 | 338010 | 338210 |
| 2 | 338020 | 338220 | |
| 3 | 338030 | 338230 | |
| 4 | 338040 | 338240 | |
| 5 | 338050 | 338250 | |
| 6 | 338060 | 338260 | |

* 建议螺栓拧紧扭矩: 35Ncm

钛基底

多功能基台临时基台



377010

* 建议螺栓拧紧扭矩: 20Ncm


多功能基台临时基台



377010

* 建议螺栓拧紧扭矩: 20Ncm


SmartBox套装 378100



易扣钛帽
内含临时衬垫



胶水垫片
378070v20




2.7kg
378020v4



1.8kg
378030v4



1.2kg
378040v4



0.6kg
378050v4



易扣阴性衬垫

冠螺丝



347010


* 建议螺栓拧紧扭矩: 20Ncm




347010

* 建议螺栓拧紧扭矩: 20Ncm

修复辅助配件



六角起子
353010 (短)
353020 (长)
353030 (短)
353040 (长)




多功能基台扳手
354050



多功能基台安装杆
353060



六角起子
353010 (短)
353020 (长)
353030 (短)
353040 (长)



多功能基台扳手
354050



多功能基台安装杆
353050



基台扳手
358010

多功能基台抛光保护帽



357020

多功能基台抛光保护帽



357020

衬垫安装工具



358020